



Big Blue® Transgenic Rodent Mutation Assay

BioReliance's Big Blue® Transgenic Rodent Mutation (TRM) assay utilizes a novel transgenic animal model bred for specific use in an assay defined by OECD Test Guideline 488, "Transgenic Rodent Somatic and Germ Cell Gene Mutation Assays." This assay allows the measurement of mutations in any tissue including germ cells and has been re-qualified for commercial use.

What is Big Blue®?


Big Blue® mice are bred to have multiple copies of recoverable target genes integrated into their genome. Big Blue® animals are created by a microinjection of the lambda shuttle vector containing the *cII* gene into the pronucleus of fertilized eggs from either C57Bl/6 mice or Fisher 344 rats.

Why the Big Blue® Transgenic Rodent Mutation Assay?

The Big Blue® TRM assay provides a means to measure mutations in any somatic cell or germ cells. Most other in vivo assays measure DNA damage, but DNA damage does not always result in mutations. The Big Blue® TRM assay thus provides the advantage of specifically measuring gene mutations in vivo and is valuable as a second tier assay to investigate positive results from in vitro gene tox assays. TRM has been included in various regulatory guidelines and expert recommendations.

What is the Big Blue® Transgenic Rodent Mutation Assay?


The Big Blue® Transgenic Rodent Mutation Assay (TRM) is an in vivo mutation assay currently available in mice. Gene mutation assays are recommended for regulatory submission, as they specifically measure gene mutations, not just DNA damage as other frequently used assays. BioReliance was the original validation laboratory for this assay in conjunction with the Institute of Environmental Health Sciences (NIEHS) and today has qualified the design in accordance with OECD Test Guideline 488.



The principle of the assay is that you can recover the lambda shuttle vector from genomic DNA and measure mutations in the lambda *cII* gene.

Big Blue® Assay Overview

In a typical study, animals are dosed by oral gavage for 28 days. Tissues are then collected at day 31 after a three day period allowing DNA lesions to be fixed into the DNA and expressed as stable mutations. After collection, high molecular weight DNA is extracted from tissues of interest and purified, the shuttle vector recovered, packaged into phage, infected onto *E. coli*, and plated for plaque formation at two temperatures. *cII* mutant phage form plaques only at 24°C while all phage form plaques at 37°C. Mutation frequency is denoted by dividing the number of mutant plaques over total number of plaques. Bioavailability (TK) measurements can be added to demonstrate exposure. OECD Test Guideline 488 allows for differing designs, so please consult a BioReliance representative to design a study that will best fit your needs and the regulatory guidelines.



Fully Qualified, Commercial Assay

Co-developed and validated by BioReliance with NIEHS in the 1990s. The assay has been re-qualified according to applicable regulatory guidelines

Approved Guideline and Regulatory Acceptance

- OECD Test Guideline 488 has been approved
- Recommended by ECHA for REACH registration
- Indicated by ICH S2(R1) and other guidelines as a follow-up to positive in vitro mutation assays

Why BioReliance and Big Blue®?

- Over 20 years experience with TRM assay by same Study Director
- Comprehensive understanding and assays for follow-on to positive test systems
- Rapid study initiation from time of inquiry
- High efficiency of mutation detection can potentially reduce the number of assays required

BioReliance Toxicology Services

Big Blue® Transgenic Rodent Mutation Assay


Re-qualification Data

Mutant Frequency Data

Representative Animal and Group (Liver)

Fold Increase

	Animal#	Phage Screened	# Mutants	Mutant Frequency (X 10 ⁻⁶)	# Packages	Group Mean and SD (X 10 ⁻⁶)
Vehicle Control	4101	978400	34	34.8	2	33.1 ± 5.4*
	4102	1750000	46	26.3	2	
	4103	916700	27	29.5	2	
	4104	1036700	36	34.7	2	
	4105	1351700	58	42.9	2	
	4106	1835100	66	36.0	3	
	4107	2311700	78	33.7	2	
	4108	1816700	49	27.0	3	
ENU Treatment	4117	1468400	156	106	2	182.0 ± 59.6*
	4118	1300000	191	147	2	
	4119	628400	167	266	2	
	4120	1185000	216	182	2	
	4121	1608300	214	133	2	
ENU Treatment	4127	658400	136	207	2	210.1 ± 40.5*
	4128	563400	156	277	2	
	4129	1136700	208	183	2	
	4130	1396700	242	173	2	
	4131	1076600	227	211	2	
ENU Treatment	4137	1398300	270	193	2	188.5 ± 46.3*
	4138	1465000	254	173	2	
	4139	830000	167	201	3	
	4140	648400	163	251	3	
	4141	1593300	197	124	2	



*Mutant frequency was homogeneous within and between groups and ENU induced statistically significant increase in *cII* mutants relative to vehicle control ($p < 0.0001$)

References

- Agilent Technologies, 2009a. λ Select-*cII* Mutation Detection System for Big Blue Rodents, Instruction Manual. Agilent Document 720120, Revision A. Santa Clara, CA. <http://www.chem-agilent.com/pdf/strata/720120.pdf>.
- Douglas, G. R., Gingerich, J. D., Gossen, J. A. & Bartlett, S. A. (1994) Sequence spectra of spontaneous lacZ gene mutations in transgenic mouse somatic and germline tissues. *Mutagenesis* 9, 451-458.
- Kohler, S.W., Provost, G.S., Fieck, A., Kretz, P.L., Bullock, W.O., Putman, D.L., Sorge, J.A. and Short, J.M. (1991) Spectra of spontaneous and mutagen-induced mutations in the lacI gene in transgenic mice. *Proc. Natl Acad. Sci. USA*, 88, 7958-7962.
- ECHA-13-R-01-EN, 2013, ECHA/2011/217, 2012
- OECD (Organisation of Economic Co-operation and Development) Guideline for Testing of Chemicals. 2011. Test Guideline 488, - Transgenic Rodent Somatic and Germ Cell Gene Mutation Assays, adopted July 28, 2011, OECD, Paris, France, <http://browse.oecdbookshop.org/oecd/pdfs/free/9748801e.pdf>.
- Piegorsch, W, W, et al. (1995) Study design and sample sizes for lacI transgenic mouse mutation assay. *Environ. and Molecular Mutagenesis*, 25: 231-245.
- Putman, D.L., Ritter A.P., Carr, G.J., Young, R.R. (1997) Evaluation of spontaneous and chemical-induced lacI mutations in germ cells from lambda/lacI transgenic mice. *Mutation Research*, 388, 129-136.
- Russell, W. L., P. R. Kelly, P. R. Hunsicker, J. W. Bangham, S. C. Maddux, and E. L. Phipps. 1979. Specific locus test shows ethylnitrosourea to be the most potent mutagen in mouse. *Proc. Natl. Acad. Sci. USA* 76:5918-5922.
- Yauk, C.L., Dubrova, Y.E., Grant, G.R., and Jeffreys, A.J. 2002. A novel single molecule analysis of spontaneous and radiation-induced mutation at a mouse tandem repeat locus, *Mutat. Res.* 500: 147-156.

www.bioreliance.com
Toll Free: 800 553 5372 Tel: 301 738 1000
Email: toxicology@bioreliance.com

 **BioReliance**[®]
by **SAFC**