

Testing Services for Veterinary Vaccines

Contents

Introduction	1
General requirements	1
Viral extraneous agents	1
Non-viral extraneous agents	1
Custom testing	1
Cell Banks	4
MCS and WCS direct testing	4
MCS and WCS extract testing	5
Mammalian Vaccines	6
Neutralizing the virus	6
<i>In vitro</i> assay	6
Avian Vaccines	7
Neutralizing the virus	8
Test for extraneous viruses using embryonated hens' eggs (Ph Eur, CVMP)	8
<i>In vitro</i> assays:	9
Test in chick kidney cells (Ph Eur, CVMP)	9
Test for <i>Avian leucosis viruses</i> (Ph Eur, CVMP)	9
Test for <i>Reticuloendotheliosis virus</i> (REV) (CVMP)	9
Test for <i>Chick anaemia virus</i> (CAA) (CVMP)	9
Test for extraneous agents using chicks (Ph Eur, CVMP)	9
Fish Vaccines	11
Vaccines Based on Insect Cell Lines/Viruses	11
Other Substances of Animal Origin	12
Preparation of test substances	12
<i>In vitro</i> assay	12
Bovine serum	12
Endpoint Tests	13
Staining for cytopathic effects	13
Haemadsorption assay	13
Immunofluorescence assay	13
Other tests	13
References	13

Introduction

General Requirements

Of paramount importance in ensuring the safety of live and inactivated veterinary vaccines is their characterization and demonstration of freedom from extraneous agents. Biological substrates, ingredients and products for veterinary vaccines produced *in vivo* and *in vitro* must be investigated for the presence of potential contaminants including; viruses, bacteria, fungi and mycoplasma in addition to stability and identity where appropriate. Both the EU and USA provide guidelines and specific recommendations on extraneous agent testing of veterinary medicinal products through the Committee for Medicinal Products for Veterinary Use (CVMP)¹⁻³, the European Pharmacopoeia (Ph Eur)⁴ and the Code of Federal Regulations (Title 9; Animals and Animal Products) (9CFR)⁵. A comprehensive testing strategy is therefore required in order to address the issues of product origin, target species and market place of the test material which could encompass:

- Master, working and cell seeds at the highest passage level
- Master virus seeds
- Final product lots
- Materials of animal origin (e.g. trypsin, serum, plasma etc)

Viral Extraneous Agents

For the detection of potential viral contaminants, suitably prepared test materials are propagated in culture to allow amplification of viral contaminants using cell lines sensitive to the viruses of concern. Depending on the test material, these cell lines will include some or all of the following:

- Primary cells of the source species (Ph Eur , CVMP)
- Cells sensitive to viruses pathogenic for the species for which the vaccine is intended (Ph Eur , CVMP, US 9CFR)
- Cells sensitive to Pestiviruses (Ph Eur , CVMP)
- Cells of the species of the cell line in which the vaccine is produced (US 9CFR)
- African green monkey kidney cells (Vero) (US 9CFR)

The choice of cells for each of the above categories is based on the sensitivity of a cell line to potential viral contaminants of the test material. The potential viral contaminants detailed in CVMP guidelines (**Table 1**) may include one or more of the following species; bovine, ovine, caprine, porcine, equine, feline, canine, rabbit, rodent and African green monkey (Vero). Potential contaminants specific to avians, fish and insects are detailed on pages 7–11. Cultures are maintained for a minimum specified period with regular subculture and observation for cytopathic effects (cpe) and morphological change. During and/or following the culture period, endpoint tests for viral contaminants are performed on the detector cells. These include general virus screening tests such as cytological staining and haemadsorption assay (HA), in addition to specific virus screening tests such as immunofluorescence assay (IFA), enzyme linked immunosorbent assay (ELISA) and polymerase chain reaction (PCR).

Non-Viral Extraneous Agents

For the detection of mycoplasma, sterility and any other non-viral species, specific testing is performed as required. The testing required is dependant on the stage of the cell seed under test, the species of origin of the cell line or virus and the target species of the product (**Table 1**).

Custom Testing

In the veterinary industry it is a common requirement that a product is licensed for both the European and US market. Accordingly all biological substrates and the ingredients used in their manufacture should be demonstrated to be free from extraneous agents using tests that satisfy distinct requirements prescribed in European and US guidelines and regulations. At BioReliance, using our expertise and comprehensive range of standards and controls, we are able to provide testing programmes tailored to your individual requirements for a broad spectrum of veterinary vaccine substrates and products.

Table 1. Extraneous agents of concern to mammalian veterinary vaccines (CVMP)

Species	Viral agent	Non-viral agent
Bovine	<i>Adenovirus</i> subgroups 1 and 2, <i>Akabane virus</i> , <i>Aujeszky's disease virus</i> , <i>Bluetongue virus</i> , <i>Epizootic hemorrhagic disease virus</i> , <i>Bovine coronavirus</i> , <i>Bovine ephemeral fever virus</i> , <i>Bovine herpesviruses</i> , <i>Bovine leukaemia virus</i> , <i>Bovine papillomaviruses</i> , <i>Bovine parvovirus</i> , <i>Bovine papular stomatitis virus</i> , <i>Pseudocowpoxvirus</i> , <i>Bovine respiratory syncytial virus</i> , <i>Bovine rotavirus</i> , <i>Bovine viral diarrhoea virus</i> , <i>Cowpox</i> , <i>Vaccinia virus</i> , <i>Foot and mouth disease virus</i> , <i>Lumpy skin disease virus</i> , <i>Malignant catarrhal fever</i> (African and European form), <i>Parainfluenza 3 virus</i> , <i>Rabies virus</i> , <i>Rift Valley fever virus</i> , <i>Rinderpest virus</i> , <i>Vesicular stomatitis virus</i>	<i>Brucella abortus</i> , <i>Leptospira</i> spp, <i>Mycobacterium tuberculosis</i> and <i>paratuberculosis</i> , <i>Mycoplasma</i> sp, <i>Salmonella</i> sp
Ovine	<i>Adenovirus</i> subgroups 1 and 2, <i>Akabane virus</i> , <i>Aujeszky's disease virus</i> , <i>Bluetongue virus</i> , <i>Epizootic hemorrhagic disease virus</i> , <i>Bovine herpesviruses</i> , <i>Bovine leukaemia virus</i> , <i>Bovine papillomaviruses</i> , <i>Bovine viral diarrhoea virus</i> , <i>Cowpox</i> , <i>Vaccinia virus</i> , <i>Foot and mouth disease virus</i> , <i>Parainfluenza 3 virus</i> , <i>Rift Valley fever virus</i> , <i>Border disease virus</i> , <i>Borna disease virus</i> , <i>Louping ill virus</i> , <i>Nairobi sheep disease virus</i> , <i>Ross River virus</i> , <i>Scrapie</i> , <i>Orf virus</i> , <i>Peste des petits ruminants</i> .	<i>Leptospira</i> spp, <i>Mycobacterium tuberculosis</i> and <i>paratuberculosis</i> , <i>Mycoplasma</i> sp, <i>Salmonella</i> sp, <i>Chlamydia ovis</i> , <i>Brucella melitensis</i>
Caprine	<i>Adenovirus</i> subgroups 1 and 2, <i>Akabane virus</i> , <i>Aujeszky's disease virus</i> , <i>Bluetongue virus</i> , <i>Epizootic hemorrhagic disease virus</i> , <i>Bovine herpesviruses</i> , <i>Cowpox</i> , <i>Vaccinia virus</i> , <i>Foot and mouth disease virus</i> , <i>Parainfluenza 3 virus</i> , <i>Rift Valley fever virus</i> , <i>Caprine arthritis encephalitis virus</i> , <i>Orf virus</i> , <i>Maedi visna virus</i> , <i>Peste des petits ruminants</i> .	<i>Leptospira</i> spp, <i>Mycobacterium tuberculosis</i> and <i>paratuberculosis</i> , <i>Mycoplasma</i> sp, <i>Salmonella</i> sp, <i>Brucella melitensis</i>
Porcine	<i>African swine fever virus</i> , <i>Aujesky's disease virus</i> , <i>Bovine viral diarrhoea virus</i> , <i>Classical swine fever virus</i> , <i>Encephalomyocarditis virus</i> , <i>Foot and mouth disease virus</i> , <i>Haemagglutinating encephalomyelitis virus</i> , <i>Transmissible gastroenteritis virus</i> and <i>Porcine respiratory coronavirus</i> , <i>Porcine adenoviruses</i> , <i>Porcine cytomegalovirus</i> , <i>Porcine enteroviruses</i> , <i>Porcine influenza virus</i> , <i>Porcine parvovirus</i> , <i>Porcine respiratory and reproductive syndrome virus</i> , <i>Porcine vesicular exanthema virus</i> , <i>Rabies virus</i> , <i>Vesicular stomatitis virus</i>	<i>Brucella suis</i> , <i>Mycoplasma hyopneumoniae</i> , <i>Mycoplasma hyorhinis</i>

Species	Viral agent	Non-viral agent
Equine	<i>African horse sickness virus, Borna disease virus, Equine arteritis virus, Equine encephalomyelitis viruses, Equine herpesviruses, Equine infectious anaemia virus, Equine influenza virus, Japanese encephalitis virus, Rabies virus, Vesicular stomatitis virus</i>	
Feline	<i>Aujeszky's disease virus, Cowpoxvirus, Feline calicivirus, Feline herpesvirus 1, Feline leukaemia virus/Feline sarcoma virus, Feline panleukopenia virus, Feline syncytia forming virus, Rabies virus</i>	<i>Chlamydia psittaci</i>
Canine	<i>Aujesky's disease virus, Canine adenoviruses 1 and 2, Canine coronavirus, Canine distemper virus, Canine herpesvirus, Canine parvovirus, Parainfluenza 2 virus, Rabies virus</i>	<i>Brucella cannis</i>
Rabbit	<i>Arenavirus, Aujeszky's disease virus, Encephalomyocarditis virus, Myxoma virus, Shope fibroma virus, Rabbit hemorrhagic disease virus, Rabies virus</i>	
Rodent	<i>Arenavirus, Encephalomyocarditis virus, Rabies virus</i>	
African green monkey (Vero)	<i>Bovine viral diarrhoea virus, Endogenous retroviruses, Reoviruses, SV40 virus, SV5 virus</i>	<i>Mycoplasma spp,</i>

Cell Banks

Master working cell seeds (MCS), working cell seeds (WCS) and WCS at the highest passage level for the production of vaccines for veterinary use must be characterized to preclude extraneous agents and adverse properties.

Testing includes microscopy, sterility, mycoplasma, viruses, species identity, karyology and tumorigenicity (Table 2). Tests are carried out for freedom from contaminating viruses directly on cultures of MCS and WCS and also by inoculating MCS and WCS extracts onto suitable detector cell lines.

MCS and WCS direct testing

Monolayers of MCS and WCS under test should be at least 70cm² (Ph Eur) or 75cm² (US 9CFR) and cultured using conditions similar to those used for preparation of the vaccine for at least 21 days (US 9CFR) or 28 days (Ph Eur) (Figure 1). Subcultures are typically every 7 days throughout the cultivation period with regular observations for evidence of cytopathic agents. At the end of the cultivation period monolayers of specified surface


area are tested for cytopathic agents by staining, haemadsorbing agents by haemadsorption assay and specified viruses by IFA or similar specific tests (see Endpoint Tests, p. 13).

Table 2. Stages of cell culture at which testing is performed

Testing required	MCS	WCS	WCS at highest passage level
General microscopy	+	+	+
Bacteria/fungi	+	+	-
Mycoplasma	+	+	-
Viruses	+	+	-
Identification of species	+	-	+ / - ¹
Karyology ¹	+	-	+
Tumorigenicity ¹	+	-	-

¹ = not required for primary cells

Figure 1: MCS and WCS Direct Testing




MCS and WCS extract testing

At the conclusion of the MCS and WCS direct test, monolayers of at least 140cm² (Ph Eur) or 75 cm² (US 9CFR) of the cells under test are freeze/thawed at least 3 times and clarified by centrifugation. Aliquots of the resultant extract are then inoculated onto monolayers of at least 75cm² (US 9CFR) of suitable detector cells. The inoculated cells are cultured for at least 14 days with at least one subculture (US 9CFR) or by preparation of freeze/thawed extracts with inoculation onto

fresh monolayers of corresponding cell type after at least 7 days (Ph Eur) (**Figure 2**). The inoculated cultures are observed regularly for evidence of cytopathic agents. After at least 14 days post inoculation (p.i.) monolayers of specified surface area are tested for cytopathic agents by staining, haemadsorbing agents by haemadsorption assay and specified viruses by IFA or similar specific tests (see Endpoint Tests, p. 13).

Figure 2: MCS and WCS Extract Testing



Mammalian Vaccines

Master seed viruses (MSV) for mammalian veterinary vaccines may be tested to preclude extraneous agents. This testing includes propagation, identity, sterility, mycoplasma and extraneous viruses. To allow testing of the virus seed or final product for extraneous viruses in the appropriate cell culture system (*in vitro* assay), the virus may require to be neutralized. Neutralization is performed to prevent infection and propagation of the virus seed in the detector cells.

Neutralizing the virus


Polyclonal or monoclonal antibody preparations containing neutralizing antibodies to the virus seed must be prepared using an antigen distinct from the virus seed isolate. The neutralizing antiserum must be shown to be free from antibodies to potential contaminants of the virus seed and also from any

non-specific inhibitory effects on the ability of contaminating viruses to infect and propagate within the cell culture. The antiserum should be used in a minimal volume to neutralize if possible at least the virus content of 10 doses of vaccine per ml. Alternative methods to neutralize or remove the test virus may be used in the absence of a suitable neutralizing antiserum.

In vitro assay

Virus (neutralized if required) is inoculated onto monolayers of suitable detector cells of at least 70cm² (Ph Eur) or 75cm² (US 9CFR) (**Figure 3**). The inoculated cells are cultured for at least 28 days (Ph Eur) or 14 days (US 9CFR). The cells are subcultured by preparation of freeze/thawed extracts with inoculation onto fresh monolayers of corresponding cell type (Ph Eur) or by at least one subculture (US 9CFR). The inoculated cultures

Figure 3: Master Seed Testing



are observed regularly for evidence of cytopathic agents. After at least 28 (Ph Eur) or 14 (US 9CFR) days p.i. monolayers of specified surface area are tested for cytopathic agents by staining, haemadsorbing agents by haemadsorption assay and specified viruses by IFA or similar specific tests (See Endpoint Tests p. 13).

Avian Vaccines


Master seed viruses (MSV) and cell seeds for avian veterinary vaccines may be tested to preclude extraneous agents. This testing includes that detailed for mammalian vaccines with some additional specifications: any chickens, embryos and tissue cultures used in production of avian vaccines must be derived from specific pathogen free (SPF) chicken flocks.

For the detection of avian derived contaminants, *in vivo* and *in ovo* assays are utilized in addition to *in vitro* assays using suitably sensitive avian detector cell lines (**Figure 4**). A full range of assays to detect avian viral agents in accordance with the CVMP and Ph Eur would include the following tests:

- Extraneous viruses using fertilized eggs
- Extraneous viruses using cell cultures
- *Avian leucosis viruses*
- *Avian reticuloendotheliosis virus*
- Extraneous agents using chicks
- *Chick anaemia virus*

Endpoints are based on observations for death and malformation (*in vivo* or *in ovo* tests) or cpe and HA (*in vitro* tests) and specific tests such as IFA and ELISA performed on material generated from the tests.

Figure 4: Test for extraneous viruses in avian vaccines



Neutralizing the virus

Polyclonal or monoclonal antibody is produced as outlined for mammalian vaccines. The neutralizing antiserum must be shown to be free from antibodies against and free from inhibitory effects on the agents listed in **Table 3** by suitably sensitive methods. Monospecific antisera for virus neutralization is not required to be tested for antibodies against any of these viruses if it can be shown that the immunizing antigen could not have been contaminated with antigens derived from that virus and if the virus is not known to infect the species of origin of the serum. Also it is not necessary to retest sera obtained from birds from SPF chicken flocks. Batches of serum must not be prepared from any passage level derived from the virus isolate used to prepare the master seed lot or from an isolate cultured in the same cell line. The antiserum should be used minimally to neutralize (if possible) at least the virus content of 10 doses of vaccine per 0.1ml (*in vitro* tests) or 0.2ml (*in vivo* tests). Alternative methods to neutralize or remove the test virus may be used in the absence of a suitable neutralizing antiserum.

Test for extraneous agents using embryonated hens' eggs

Virus (neutralized if required) containing at least 10 doses of vaccine in 0.2ml is inoculated into 3 groups of 10 embryonated hens' eggs as follows: group 1- the allantoic cavity, group 2 – the chorio-allantoic membrane (CAM) and group 3 – the yolk sac. The eggs are candled daily for 7 days (group 1 and 2) or 12 days (group 3). All embryos that die after at least 24 hours or survive are examined macroscopically for abnormalities. The CAM of these eggs is also examined and the allantoic fluids are tested for haemagglutinating agents. A further embryo passage is performed by pooling separately the material from live, dead and abnormal embryos. Each pool is inoculated into 10 eggs for each route as detailed previously: CAM material is inoculated into CAMs, allantoic fluids into the allantoic cavity and embryo material into the yolk sac. The inoculated eggs are candled daily as detailed previously and at the end of the test period the eggs are examined and tested as for the initial period.

Table 3. Antibody specifications for antiserum used to neutralize avian vaccines (Ph Eur, CVMP)

<i>Avian adenovirus</i>	<i>Chick anaemia virus</i>
<i>Avian encephalomyelitis</i>	<i>Duck enteritis virus</i>
<i>Avian infectious bronchitis viruses</i>	<i>Duck hepatitis virus type 1 and type 2</i>
<i>Avian infectious bursal disease virus types 1 and 2</i>	<i>Egg drop syndrome virus</i>
<i>Avian infectious haemorrhagic enteritis virus</i>	<i>Fowl pox virus</i>
<i>Avian infectious laryngotracheitis virus</i>	<i>Influenza viruses</i>
<i>Avian leucosis viruses</i>	<i>Marek's disease virus</i>
<i>Avian nephritis virus</i>	<i>Turkey herpesvirus</i>
<i>Avian paramyxoviruses 1-9</i>	<i>Turkey rhinotracheitis virus</i>
<i>Avian orthoreoviruses</i>	<i>Newcastle disease virus</i>
<i>Avian reticuloendotheliosis virus</i>	

***In vitro* assays**

Test in chick kidney cells (Ph Eur, CVMP)

Virus (neutralized if required) is inoculated onto 5 replicate monolayers of at least 25cm² of chick kidney cells. The inoculated cells are cultured for at least 21 days with subculture at 4 to 7 day intervals. Each subculture is performed with pooled cells and fluids from all 5 monolayers that have undergone one freeze/thaw cycle. The extract is then inoculated onto fresh monolayers of chick kidney cells. The inoculated cultures are observed regularly for evidence of cytopathic agents. At the end of the culture period monolayers of specified surface area are tested for cytopathic agents by staining, haemadsorbing agents by haemadsorption assay and for haemagglutinating agents by haemagglutination assay.

Test for *Avian leucosis viruses* (Ph Eur, CVMP)

Virus (neutralized if required) is inoculated onto 5 replicate monolayers of at least 50cm² of primary or secondary chick embryo fibroblasts that are known to be susceptible to subgroups A, B and J of *Avian leucosis viruses* (support the growth of exogenous but not endogenous *Avian leucosis viruses*). The inoculated cells are cultured for at least 9 days with subculture at 3 to 4 day intervals. Cells are retained from each subculture and at the end of the culture period are tested for group specific *Avian leucosis* antigen by ELISA assay.

Test for *Reticuloendotheliosus virus* (REV) (CVMP)

Virus (neutralized if required) is inoculated onto 5 replicate monolayers of at least 25cm² of primary or secondary chick or duck embryo fibroblasts. The inoculated cells are cultured for at least 10 days with subculture twice at 3 to 4 day intervals. At the end of the culture period monolayers of specified surface area are tested for REV by IFA.

Test for *Chick anaemia virus*

DNA is extracted from an appropriate volume of virus and tested by polymerase chain reaction (PCR) to detect *Chick anaemia virus*. The PCR assay allows the detection of a nucleic acid target molecule and uses two target specific oligonucleotide primers that flank a target DNA sequence. The internal sequence is amplified by repeated cycles of heat denaturation of the DNA template, annealing the primers to their complementary sequences on each DNA strand and extension of the annealed primers with a thermostable Taq DNA polymerase. The amplified sequences are identified by hybridization with an oligonucleotide probe specific for the target DNA.

Test for extraneous agents using chicks (Ph Eur, CVMP)

At least 10 chicks that are 2 weeks old (older birds may be used if the seed virus is pathogenic for birds of this age) are inoculated with virus (neutralized if required). Virus is inoculated at 100 doses by the intramuscular route and 10 doses by eye-drop. Repeat inoculations are performed 2 weeks later. The chicks are observed for a period of 5 weeks. Serum is collected from each chick at the end of the test period and is tested for antibodies to extraneous agents using specific tests listed in **Table 4** and also those listed in **Tables 5, 6 and 7** if required.

Table 4. Standard antibody tests for extraneous agents using chicks (Ph Eur, CVMP)

Agent	Detection test
<i>Avian adenoviruses group 1</i>	SN, EIA, AGP
<i>Avian encephalomyelitis virus</i>	AGP, EIA
<i>Avian infectious bronchitis virus</i>	EIA, HI
<i>Avian infectious laryngotracheitis virus</i>	SN, EIA, IS
<i>Avian leucosis viruses</i>	SN, EIA
<i>Avian nephritis virus</i>	IS
<i>Avian orthoreoviruses</i>	IS, EIA
<i>Avian reticuloendotheliosis virus</i>	AGP, IS, EIA
<i>Chick anaemia virus</i>	IS, EIA, SN
<i>Egg drop syndrome virus</i>	HI, EIA
<i>Avian infectious bursal disease virus</i>	Serotype 1 - AGP, EIA, SN, serotype 2 – SN, AGP, EIA, HI
<i>Influenza A virus</i>	AGP, EIA, HI
<i>Marek's disease virus</i>	AGP
<i>Newcastle disease virus</i>	HI, EIA
<i>Turkey rhinotracheitis virus</i>	EIA
<i>Salmonella pullorum</i>	Agg

Agg: agglutination

IS: immunostaining

AGP: agar gel precipitation

HI: haemagglutination inhibition

EIA: enzyme immunoassay

SN: serum neutralization

If the seed virus is of turkey origin or was propagated in turkey substrates, tests for the following agents are also performed:

Table 5. Additional antibody tests for turkey extraneous agents

Agent	Detection test
<i>Chlamydia spp.</i>	EIA (CFT or AGP) ¹
<i>Avian infectious haemorrhagic enteritis virus</i>	AGP (FAT or VN) ¹
<i>Avian paramyxovirus 3</i>	HI (ELISA) ¹
<i>Avian infectious bursal disease virus type 2</i>	SN
<i>Turkey lympho-proliferative disease virus</i>	20 turkey poults are inoculated by the intraperitoneal route. Sections of spleen and thymus are taken from 10 of the poults 2 weeks after inoculation and observed for macroscopic and microscopic lesions. The remaining poults are observed for at least 40 days.

¹Tests in parentheses recommended by CVMP.

CFT: complement fixation test VN: virus neutralization

FAT: fluorescence antibody test

If the seed virus is of duck origin or was propagated in duck substrates, tests for the following agents are also performed:

Table 6. Additional antibody tests for duck extraneous agents

Agent	Detection test
<i>Chlamydia</i> spp.	EIA (CFT or AGP) ¹
Duck and goose parvovirus	SN, EIA
Duck enteritis virus	SN (ELISA) ¹
Duck hepatitis virus type 1 (and 2) ¹	SN (FAT) ¹

¹ Tests in parentheses recommended by CVMP.

If the seed virus is of goose origin or was propagated in goose substrates, tests for the following agents are also performed:

Table 7. Additional antibody tests for goose extraneous agents

Agent	Detection test
Duck and goose parvovirus	SN, EIA
Duck enteritis virus	SN
Goose haemorrhagic polyomavirus	At least 10 doses of virus (neutralized if required) is inoculated subcutaneously to each of 10 susceptible goslings. The gosling are observed for at least 28 days.

Fish Vaccines

Master seed viruses (MSV) and cell seeds for fish veterinary vaccines may be tested to preclude extraneous agents. This testing includes those detailed for mammalian vaccines with some additional specifications as detailed in **Table 8**.

Table 8. Specific agents for fish vaccines

Viral agents	Protozoal agents	Bacterial agents
Viral haemorrhagic septicaemia virus (VHSV)	<i>Myxosma cerebralis</i>	<i>Yersinia ruckeri</i>
Infectious haematopoietic necrosis virus (IHNV)		<i>Vibrio anguillarum</i>
Spring viremia of carp virus (SVCV)		<i>Aeromonas salmonicinarum</i>
Infectious pancreatic necrosis virus (IPNV)		

Vaccines based on Insect Cell Lines/Viruses

Master seed viruses (MSV) and cell seeds of insect origin may be tested to preclude extraneous agents. This testing includes those detailed for mammalian vaccines with some additional testing for potential contaminants of insect origin:

- Appropriate detector cell lines including one or more of insect origin and one or more capable of detecting extraneous viral agents of insect origin
- *In vivo* assay in suckling mice to detect virus, including *Arbovirus*, that may not cause cpe or other discernable effects *in vitro*
- Electron microscopy assay for the detection and quantitation of viral particles
- F-PERT assay for the detection of retroviruses on cell supernatant from exponentially growing cells
- Test for spiroplasma and mycoplasma

Other Substances of Animal Origin

Substances of animal origin such as serum, trypsin and plasma may be used during the manufacture of veterinary immunological products as ingredients of culture medium etc or added constituents of vaccines or diluents. These materials must be tested for freedom from extraneous agents once they have been prepared in an appropriate manner for inoculation onto detector cell cultures (*in vitro* assay).

Preparation of test substances (including bovine serum)

For the preparation of test substances solids are dissolved or resuspended in a suitable medium to create a solution or suspension containing at least 300g/L of the test material (Ph Eur) or 15% of the ingredient is used in the growth medium (US 9CFR).

In vitro assay for test substances (excluding bovine serum)

The test material (prepared appropriately), is inoculated onto monolayers of at least 75cm² of suitable detector cells. The inoculated cells are cultured for a total of at least 21 days (Ph Eur

and US 9CFR) with regular observations for cytopathic agents. At the end of each 7 day period a proportion of the original cultures are fixed, stained and examined for cpe, tested by HA and tested for specific agents by appropriate serodiagnostic techniques such as IFA. The remaining cells are subcultured (Ph Eur).

Alternatively cells are subcultured at least 2 times over the 21 day period at the end of which monolayers of specified surface area are tested for cytopathic agents by staining and tested by HA and IFA or similar for specified viruses (US 9CFR).

Testing of bovine serum used in the production of immunological veterinary and medicinal products

General and specific tests should be performed prior to any inactivation treatment and further tests performed post inactivation. The general and specific tests must be capable of detecting the viruses listed in **Table 9**. Testing is performed as per standard tests for mammalian vaccines and in accordance with the latest CVMP guidelines⁶.

Table 9. Testing required for bovine serum

<i>Bovine adenovirus</i>
<i>Bovine viral diarrhoea virus (BVDV)</i>
<ul style="list-style-type: none"> To test for BVDV the test serum is incorporated in culture medium used for cultivation of bovine cells sensitive to <i>Pestiviruses</i>. After at least 4 subcultures the cells are tested by IFA Serum should also be shown to be free from BVDV antibodies or that the level present does not interfere with detection of a low titre of BVDV Serum should be tested for BVDV post inactivation (this is not required if BVDV was not detected prior to inactivation) Comparative titration results – post inactivation treatment, a reference strain of BVDV should be titrated in sensitive cells that have been grown in the presence of the test serum for at least 3 subcultures.
<i>Parvovirus</i>
<i>Bovine respiratory syncytial virus</i>
<i>Reovirus</i>
<i>Parainfluenza 3</i>
<i>Infectious bronchitis virus (BHV-1)</i>
<i>Bluetongue virus</i>

Endpoint Tests

Staining for cytopathic effects

Cells infected with cytopathic viruses display cpe or morphological changes. Cells exhibiting these changes can be fixed and stained with an appropriate cytological stain. Stained cell cultures are then observed for inclusion bodies, giant cells or other abnormalities or lesions indicative of changes in cellular morphology which may be attributable to a viral contaminant (Figure 5).

Haemadsorption assay

Some members of the virus groups of orthomyxo-, toga-, parvo- and paramyxoviruses may not always produce an obvious cpe in susceptible cell cultures. These viruses may be detected by a haemadsorption assay that utilises their haemagglutinin expression. A suspension of red blood cells including chicken, guinea pig, human and any other species of concern is incubated with the test culture. After incubation at various temperatures the cultures are washed and observed microscopically to detect haemadsorption (Figure 6).

Immunofluorescence assay

Monoclonal or polyclonal antibodies specific for particular antigens are utilized to detect specific viruses in fixed cell cultures. The antibodies are labelled with a fluorescent dye and when in contact with their specific antigen the reaction site can be observed using a fluorescent microscope. If viral antigen is not present the antibodies are removed during the washing steps. This test is ideal for observing non-cytopathic and non-haemadsorbing viruses. The specific viruses to be tested for by IFA are determined by the species of origin of the cell line or virus and the intended species for the product (Figure 7).

Other tests

Other specific tests may be employed to determine the presence of certain groups of, or specific viruses. These may include physical tests such as electron microscopy (EM), *in vivo* tests in eggs or animals and chemical tests such as ELISA, neutralization, nucleic acid hybridization, reverse transcriptase assay (RT) and PCR. These tests may also be used to investigate anomalous results obtained from any of the three standard endpoint tests.

Figure 5: infected cells stained with a cytological stain





Figure 6: cells displaying haemadsorption



Figure 7: cells displaying immunofluorescence



References

- ¹ The rules governing medicinal products in the European Union, volume 7B, Guidelines, Veterinary Medicinal Products, Immunologicals, quality. General Requirements for the Production and Control of Live Mammalian Bacterial and Viral Vaccines for Veterinary Use. 7Blm1a (111/3182/91-EN), March 1992.
- ² The rules governing medicinal products in the European Union, volume 7B, Guidelines, Veterinary Medicinal Products, Immunologicals, quality. General Requirements for the Production and Control of Inactivated Mammalian Bacterial and Viral Vaccines for Veterinary Use. 7Blm1b (111/3181/91-EN), March 1992.
- ³ Specific Requirements for the Production and Control of Live and Inactivated Vaccines Intended for Fish. 7Blm9a (111/3590-92) September 1994.
- ⁴ European Pharmacopoeia (2008)
- ⁵ Code of Federal Regulations Title 9, Animals and Animal Products, June 30 2008
- ⁶ Guideline on "Requirements and controls applied to bovine serum used in the production of immunological veterinary medicinal products, EMEA/CVMP/743/00-FINAL". Committee for Veterinary Medicinal Products 2002.

BioReliance Corp.

14920 Broschart Road
Rockville, Maryland 20850
Tel: 800.553.5372
Fax: 301.610.2590
Email: info@bioreliance.com

BioReliance Ltd.

Todd Campus
West of Scotland Science Park
Glasgow, Scotland G20 0XA
Tel: 44 (0) 141 946.9999
Fax: 44 (0) 141 946.0000
Email: info@bioreliance.com

BioReliance Ltd.

Innovation Park
Hillfoots Road
Stirling, Scotland FK9 4NF
Tel: 44 (0) 141 946.9999
Fax: 44 (0) 141 946.0000
Email: info@bioreliance.com

BioReliance, K.K.

c/o Sigma-Aldrich Japan K.K.
Tennoz Central Tower 4F
2-2-24 Higashi-Shinagawa
Shinagawa-ku
Tokyo 140-0002, Japan
Tel: +81 (0)3 5796 7430
Fax: +81 (0)3 5796 7435
Email: jpinfo@bioreliance.com

BioReliance Ltd.

c/o Sigma-Aldrich Chemicals Pvt Ltd
102 Alpha Building
Hiranandani Gardens
Powai, Mumbai 400076
Tel: +91 22 40872364
Fax: +91 22 25797589
Email: indiainfo@bioreliance.com

www.bioreliance.com

North America Toll Free: 800 553 5372 Tel: 301 738 1000
Europe & International Tel: +44 (0)141 946 9999 • Japan Tel: +03 5796 7430
Email: info@bioreliance.com